



Cleo II

Advanced technology, improved ergonomics and the Operator-friendly features of the Cleo II ensure superior control and precision for the dental team. The flexibility of the Cleo II range ensures that any configuration will equip a dental surgery with the ultimate, complete surgery system. The unique folding and extending chair provides unrivalled comfort and access for both Operator and Patient.

CLEO II 'Surgery System'

Air or Electric versions



CLEO II 'Flexible System'¹ (A) or (E)



Mobile Cart
Air or Electric versions



Cabinet Module
Air or Electric versions

¹ Available separately / Cleo II operating light options, see page 15

Patient-friendly, synchronised chair :

Unique folding-extending chair occupies an unusually small footprint and provides unrivalled comfort, natural support and enhanced access for both Operator and Patient.

Seamless upholstery

Twin-axis backrest with lumbar support

Twin-articulating headrest

Folding and extending legrest

Seamless, easy-clean upholstery offers **optimum hygiene levels**

Foldaway, detachable armrest

Wheelchair access is made possible by the unique combination of an inclined chair and opened armrest

A folded legrest and foldaway armrest ensures zero obstruction to Patient entry and exit

Folding and extending legrest safety feature ensures that the legrest retracts if contact is made with any object

uniquely ideal :

Cleo II

Equipped for an array of surgery procedures and working methods in complete comfort.

Height-adjustable, extendable, rotating Assistant's console

Easy-clean porcelain spittoon with cup-filler and bowl-flush

Rotating, chair-mounted cuspidor (90°)

Rotating cuspidor unit increases access for an Assistant working closely over the Patient from the 3 o'clock position

'E' version key features ~

- LED digital display
- Preset micro motor
- Torque control with Auto-Reverse
- First Priority handpiece selection

Operator's console is height-adjustable and rotates: supreme versatility and the pneumatic brake enables precise positioning for ideal working comfort.

Removable rubber mat



Control and precision are the prerequisites of the Cleo II. The Operator's console rotates smoothly on a compact, base-mount enabling an Operator easy access between the 9 and 1 o'clock positions. Hygiene-conscious at every point, the Cleo II is designed to perform in the most demanding dental surgery environments.

Versatility for **working comfort** :

From consultation through treatment, Cleo II provides supreme access and working comfort for the Operator and Assistant.




Folding legrest and detachable, foldaway armrest for perfect Patient entry and exit

Clean and preparation during consultation improves efficiency



Maximised workflow and a greater variety of operation are the key achievements of the Cleo II treatment centre. An inclined Cleo II chair is a very gentle welcome to any Patient entering the surgery and enables the most pleasant consultation position possible.



complete, hygienic surgery system :

Cleo II




Operator can work unobstructed from the front, side or behind the Patient.

CLEO II 'Surgery System'
permits multiple dental disciplines
with increased Operator comfort

Ideal for **single-** and **four-handed**
working positions

'E' version



Extending and rotating
the Assistant's console
enables very comfortable
single-handed operation

*Rotating cuspidor unit increases
access for an Assistant working
closely over the Patient from the
3 o'clock position*

Cleo II 'Flexible' treatment centres :

Mobile configurations to suit any Practice and the flexibility of Electric and Air versions to fulfil individual preferences.



CLEO II 'FLEXIBLE SYSTEM'

Tubes are easily detached

Rotating cuspidor (90°)

Foldaway, detachable armrest

Vacuum solids collector is easily accessed



Ultra-slim backrest improves Operator's working comfort

Integrated foot control stalks for chair movements including two programmable presets, Last-Position memory and Auto-Return

Cleo II chair key features ~

- Folding and extending legrest
- Twin-articulating headrest
- Foldaway, detachable armrest
- Synchronised chair movements
- Lumbar support
- Seamless upholstery

Folding and extending legrest safety feature ensures that the legrest retracts if contact is made with any object

The telescopic base is a superior, elegant design solution for easy hygiene maintenance – smooth, white surfaces, with integrated foot controls. Supremely comfortable chair movements, smooth, near-silent and gentle start/stop.

technically advantaged :

Cleo II

Super ergonomics for Dentist and Patient.

'E' version control panel includes instrument, chair, spittoon and light controls



'A' version control panel includes chair and spittoon controls



Increased range of synchronised chair movements



Mobile Cart and Cabinet Module 'A' and 'E' versions



Assistant's console with chair, spittoon and light¹ controls

Control all chair movements by hand or foot: by hand from both the Operator's and Assistant's console; by foot using the stalks built-in to the chair base. Spittoon functions and ON/OFF light control ¹ are also activated from both the Operator's and Assistant's console.



820S LIGHT – Unit-mounted or 700 SERIES LIGHT – Unit- or ceiling-mounted

Triple-axis head movement
White natural beam 5,000 Kelvin
Touchless sensor ON/OFF control
Composite Curing Safe Mode

4,200 Kelvin
Touchless sensor ON/OFF control
Composite Curing Safe Mode
Retractable Patient mirror

Operator's foot control for instruments

¹ Light controls on 'E' version only

Ceiling-mounted NDL and 700 Series lights also available, see 'Operating Lights' on pages 48-51