



Stools



ic-Tilt



ic-DF



ic-D



ic-DSSA



ic-DSAF



ic-DSA

A comprehensive range of dental stools for Operators and Assistants. Ergonomic designs to provide full support in the most demanding working conditions. Intelligent features such as tilting and rotating backrests to improve posture and simple pleasures like a soft-cushioned seat are essential for comfort during operation.

Each stool uses a long-lasting, hydraulic lift for smooth height-adjustment and a five-point base with smooth-running castors. Standard and Ultrasoft upholstery is available in a vast array of colours¹ to suit any taste and complete any surgery design.

Operators' Stools

Operators' stools feature tilting and fixed backrest options.

ic-TILT Tilting and height-adjustable backrest



Operators' stools ~ 465mm min. seat height / 590mm max. seat height

ic-DF – with ergonomic seat and (fixed) backrest



ic-D – the simpler Operators' stool with fixed, centred backrest



Ergonomic seat and backrest are constructed with dual-density cushioning to provide the best support for each part of the body



Assistants' Stools

Assistants' stools feature a footring, with rotating and fixed backrest options.

ic-DSAF – the alternative Assistants' stool with ergonomic seat and backrest



ic-DSA – with fixed, offset backrest



Soft-cushioned seat and backrest



ic-DSSA

Rotating backrest for support



Assistants' stools ~ 530mm min. seat height / 710mm max. seat height