

PERFORMER[®]



a dec[®]

The A-dec Difference

Legendary Reliability

A-dec has earned a reputation for creating equipment of legendary reliability and quality. It is a reputation earned through dedication to excellence, and commitment to innovative, yet simple engineering design. This means you can count on easy maintenance and a long product life.



Control block – Patented design promotes clean treatment water with every activation.



Lift cylinder – A-dec built and engineered with robust hydraulics and fewer parts.



A-dec LED light – Award-winning silent design incorporates unique “passive cooling” to dissipate heat without a fan.

Assurance

With the Performer, you are assured comfort and simplicity. Add to that the



peace of mind that comes with A-dec's five-year warranty. It's all part of what makes A-dec, A-dec.



Shown with Performer light,
monitor and back-mount delivery

Performer

Features Overview


Designed to give you flexibility, access and comfort, the Performer answers the need for an entry-level option, plus adds a rich feature set that supports productivity.

Shown with Performer light, cuspidor, monitor and back-mount delivery



Choices

The Performer is known for its reliability. Now with added functionality and a multitude of options, you can feel secure in the knowledge that this is a chair built for the long term.



Front- or back-mounted delivery system



Disc or Lever foot control



Standard or Deluxe touchpad



Monitor mount available



Rotating vitreous china cuspidor available



Assistant's 3-position holder with Durr HVE, saliva ejector switched holders and optional syringe




Performer light

6300 Halogen light

LED light

Left- or Right-handed, You Decide

Performer's delivery system and assistant's instrumentation is easy to switch between right- and left-handed users. The ambidextrous design offers full flexibility and efficiency for the team, and creates a more flexible setup for patient entry and exit.



Seamless Upholstery Color Choices



Order Your Samples at
a-dec.com/InspireMe



Shown with Performer light,
monitor and back-mount delivery

Quick Facts

- Left- or right-hand configuration
- Smooth, quiet hydraulic cantilever system
- Choice of Performer light, 6300 halogen light or LED light
- Choose a front- or back-mounting location for the delivery system
- Large 2-liter water bottle requires fewer refills
- Monitor mount with swivel and tilt
- User programmable touchpad
- Assistant's instrumentation SYR/HVE/SE
- Safety lift arm stop plate
- Robust cast metal construction
- 400 pound (181kg) patient capacity
- Single or two-piece floorbox accommodates placement of utilities
- Contemporary and stylish delivery system
- 300 watt power supply

Meet the Rest of the Family

Along with the Performer, A-dec offers a complete line of dental equipment solutions. If you are looking for additional premium features, we invite you to learn more about the rest of our product family, from the A-dec 300, A-dec 400 and A-dec 500 dental chairs, to delivery and support systems.



A-dec 500



A-dec 400



A-dec 300



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132 USA
Tel: 1.800.547.1883 within USA/CAN
Tel: 1.503.538.7478 outside USA/CAN
Fax: 1.503.538.0276
a-dec.com

A-dec Australia

Tel: 1.800.225.010 within AUS
Tel: +61.(0)2.8332.4000 outside AUS

A-dec China

Tel: +1.503.538.7478

A-dec United Kingdom

Tel: 0800.233285 within UK
Tel: +44.(0).24.7635.0901 outside UK

A-dec, the A-dec logo, A-dec 500, A-dec 300, Cascade, Cascade Master Series, Century Plus, Continental, Decade, ICX, ICV, Performer, Preference, Preference Collection, Preference ICC, and Radius are trademarks of A-dec Inc. and are registered in the United States and other countries. A-dec 400, A-dec 200, Preference Slimline, and reliablecreativesolutions are also trademarks of A-dec Inc. None of the trademarks or trade names in this document may be reproduced, copied, or manipulated in any manner without the express, written approval of the trademark owner.

Certain touchpad symbols are proprietary to A-dec Inc. Any use of these symbols, in whole or in part, without the express written consent of A-dec Inc., is strictly prohibited.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.